


Rabbit Anti-Human Pth1R Monoclonal Antibody For Flow Cytometry - Q03431

Item Number: CM0009955



Introduction

-conjugated rabbit monoclonal antibody targeting human PTH1R (parathyroid hormone 1 receptor), validated for flow cytometry. Suitable for GPCR signaling and bone/kidney disease research. Specific to human PTH1R.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human PTH1R mAb
Remarks/Alias	PFE; EKNS; PTHR; PTHR1
Species	Human
GeneID (human)	5745
GeneID	5745
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-188 of human PTH1R (NP_000307.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q03431
Protein Weight	66kDa
Shipping	Ice bag