


Rabbit Anti-Human Ror1 Monoclonal Antibody For Flow Cytometry - Q01973

Item Number: CM0008137



Introduction

Conjugated rabbit monoclonal antibody targeting human ROR1, validated for flow cytometry. Specific to human ROR1, it detects a 104kDa protein. Suitable for Wnt signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human ROR1 mAb
Remarks/Alias	NTRKR1; dj537F10.1
Species	Human
GeneID (Human)	4919
GeneID	4919
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-403 of Human ROR1 (NP_005003.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q01973
Protein Weight	104kDa
Shipping	Ice bag