

## Rabbit Anti-Human Siglec-8 Monoclonal Antibody For Fc - Q9Nyz4

Item Number: CM0011688



### Introduction

-conjugated rabbit monoclonal antibody targeting human Siglec-8 (SAF-2). Validated for flow cytometry, it detects the 44-54 kDa membrane lectin on eosinophils, basophils, and mast cells. Ideal for immunophenotyping and allergy research.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human Siglec-8 mAb                                                                                       |
| Remarks/Alias    | SAF2; SIGLEC-8; SIGLEC8L                                                                                             |
| Species          | Human                                                                                                                |
| GeneID           | 27181 (Human)                                                                                                        |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 17-363 of human Siglec-8 (NP_055257.2) |
| Source           | Rabbit                                                                                                               |
| Category         | Monoclonal Antibodies                                                                                                |
| Application      | FC                                                                                                                   |
| Cross Reactivity | Human                                                                                                                |
| SWISS            | Q9NYZ4                                                                                                               |
| Protein Weight   | 44kDa/46kDa/54kDa                                                                                                    |
| Shipping         | Ice bag                                                                                                              |