

Rabbit Anti-Human Siglec-8 Monoclonal Antibody For Flow Cytometry - Q9Nyz4

Item Number: CM0011530



Introduction

-conjugated rabbit monoclonal antibody targeting human Siglec-8 (Q9NYZ4), validated for flow cytometry. Specifically detects Siglec-8 on basophils, mast cells and eosinophils.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Siglec-8 mAb
Remarks/Alias	SIGLEC8, SAF2, SIGLEC-8, SIGLEC8L, sialic acid-binding Ig-like lectin 8
Species	Human
Gene ID	27181
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-363 of human Siglec-8 (NP_055257.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9NYZ4
Protein Weight	44kDa/46kDa/54kDa
Shipping	Ice bag