

Rabbit Anti-Human Siglec-8 Monoclonal Antibody For Flow Cytometry - Q9Nyz4

Item Number: CM0009340



Introduction

-conjugated rabbit monoclonal antibody targeting human Siglec-8 (SAF-2). Designed for flow cytometry to detect eosinophils, basophils, and mast cells. Recombinant immunogen ensures high specificity.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Siglec-8 mAb
Remarks/Alias	SAF2; SIGLEC-8; SIGLEC8L
Species	Human
GeneID	27181 (human)
Immunogen	Recombinant protein (a.a. 17-363 of human Siglec-8, NP_055257.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9NYZ4
Protein Weight	44kDa/46kDa/54kDa
Shipping	Ice bag