

Rabbit Anti-Human Siglec9 Monoclonal Antibody For Flow Cytometry - Q9Y336

Item Number: CM0009499

Introduction

-conjugated rabbit anti-human SIGLEC9 monoclonal antibody for flow cytometry. Detects CD329 on human neutrophils and monocytes. Recombinant immunogen; protein weight 50-52 kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human SIGLEC9 mAb
Remarks/Alias	CD329; CDw329; FOAP-9; siglec-9; OBBP-LIKE
Species	Human
Gene ID	27180 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-348 of human SIGLEC9 (NP_055256.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9Y336
Protein Weight	50kDa, 52kDa
Shipping	Ice bag