

Rabbit Anti-Human Syk Mab For Fc (Intra) - P43405

Item Number: CM0006253



Introduction

-conjugated rabbit monoclonal antibody targeting human Syk (p72-Syk). Designed for intracellular flow cytometry. Ideal for investigating B-cell receptor signaling, innate immunity, and platelet activation pathways.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Syk mAb
Remarks/Alias	IMD82; p72-Syk
Species	Human
GeneID (Human)	6850
GeneID	6850
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 251-350 of human Syk (NP_003168.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P43405
Protein Weight	72kDa/70kDa
Shipping	Ice bag