

## Rabbit Anti-Human Syk Monoclonal Antibody For Fc (Intra) - P43405

Item Number: CM0010992



### Introduction

-conjugated rabbit monoclonal antibody targeting human Tyrosine-protein kinase SYK (p72-Syk). Validated for intracellular flow cytometry in human samples. Essential for B-cell receptor signaling and innate immunity studies. UniProt P43405.

[Learn More](#)

| Parameter        | Value                                                                                               |
|------------------|-----------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human Syk mAb                                                                           |
| Remarks/Alias    | IMD82; p72-Syk                                                                                      |
| Species          | Human                                                                                               |
| Gene ID (Human)  | 6850                                                                                                |
| Gene ID          | 6850                                                                                                |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 251-350 of human Syk (NP_003168.2) |
| Source           | Rabbit                                                                                              |
| Category         | Monoclonal Antibodies                                                                               |
| Application      | FC (intra)                                                                                          |
| Cross Reactivity | Human                                                                                               |
| SWISS            | P43405                                                                                              |
| Protein Weight   | 72kDa/70kDa                                                                                         |
| Shipping         | Ice bag                                                                                             |