

Rabbit Anti-Human Syk Monoclonal Antibody For Intracellular Flow Cytometry - P43405

Item Number: CM0011130



Introduction

Rabbit monoclonal antibody against human Syk (Spleen tyrosine kinase) conjugated to . Suitable for intracellular flow cytometry (FC) on human samples. Syk is a key non-receptor kinase in B-cell receptor signaling and immune responses.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Syk mAb
Remarks/Alias	IMD82; p72-Syk
Species	Human
GeneID (Human)	6850
GeneID	6850
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 251-350 of human Syk (NP_003168.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P43405
Protein Weight	72kDa/70kDa
Shipping	Ice bag