

Rabbit Anti-Human Syndecan 4 Monoclonal Antibody For Fc - P31431

Item Number: CM0010381



Introduction

Rabbit anti-Human Syndecan 4 monoclonal antibody, FC-validated. Recognizes human Syndecan-4 (SYND4), a cell surface proteoglycan regulating exosome biogenesis. Recombinant immunogen. Ideal for flow cytometry.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Syndecan 4 mAb
Remarks/Alias	SYND4
Species	Human
GeneID	6385
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-145 of human Syndecan 4 (NP_002990.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P31431
Protein Weight	22kDa
Shipping	Ice bag