

Rabbit Anti-Human T-Bet/Tbx21 Monoclonal Antibody For Flow Cytometry - Q9UI17

Item Number: CM0011625

Introduction

-conjugated rabbit monoclonal antibody against human T-bet/Tbx21, validated for intracellular flow cytometry. Ideal for Th1 cell differentiation studies, immunophenotyping, and transcription factor analysis.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human T-bet/Tbx21 mAb
Remarks/Alias	TBX21; T-PET; T-bet; TBET; TBLYM; T-box 21
Species	Human
Gene ID (Human)	30009
Gene ID	30009
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 322-535 of human T-bet/Tbx21 (NP_037483.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	Q9UL17
Protein Weight	58kDa
Shipping	Ice bag