


Rabbit Anti-Human Tie-2/Cd202B Monoclonal Antibody For Fc - Q02763

Item Number: CM0004933



Introduction

-conjugated rabbit monoclonal antibody targeting human Tie-2/CD202b (Angiopoietin-1 receptor), validated for flow cytometry (FC). Ideal for angiogenesis, vascular biology, and endothelial cell research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Tie-2/CD202b mAb
Remarks/Alias	TIE2; VMCM; GLC3E; TIE-2; VMCM1; CD202B
Species	Human
GeneID (Human)	7010
GeneID	7010
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 23-745 of human Tie-2/CD202b (NP_000450.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q02763
Protein Weight	126kDa
Shipping	Ice bag