


Rabbit Anti-Human Tigit Monoclonal Antibody For Flow Cytometry - Q495A1

Item Number: CM0010639



Introduction

-conjugated rabbit monoclonal antibody against human TIGIT (Q495A1), validated for flow cytometry. Ideal for immune checkpoint research, T-cell and NK cell inhibition studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human TIGIT mAb
Remarks/Alias	VSIG9; VSTM3; WUCAM
Species	Human
Gene ID (Human)	201633
Immunogen	Recombinant protein or recombinant fusion protein containing a sequence corresponding to amino acids 22-141 of human TIGIT (NP_776160.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q495A1
Protein Weight	18kDa/26kDa
Shipping	Ice bag