


Rabbit Anti-Human Tim-3/Havcr2 Monoclonal Antibody For Flow Cytometry - Q8Tdq0

Item Number: CM0009771



Introduction

-conjugated rabbit monoclonal antibody targeting human TIM-3/HAVCR2, validated for flow cytometry (FC). Ideal for studying immune checkpoint regulation, T-cell exhaustion, and cancer immunotherapy research.

[Learn More](#)

Specification	Detail
Product Name	Rabbit anti-Human TIM-3/HAVCR2 mAb
Remarks/Alias	TIM3; CD366; KIM-3; SPTCL; TIMD3; Tim-3; TIMD-3; HAVcr-2
Species	Human
GeneID	84868
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8TDQ0
Protein Weight	33kDa
Shipping	Ice bag