


Rabbit Anti-Human Tim-3/Havcr2 Monoclonal Antibody For Flow Cytometry - Q8Tdq0

Item Number: CM0005164



Introduction

This -conjugated rabbit monoclonal antibody targets human TIM-3/HAVCR2 (CD366), ideal for flow cytometry analysis of immune checkpoint expression on T cells, NK cells, and dendritic cells.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human TIM-3/HAVCR2 mAb
Remarks/Alias	TIM3; CD366; KIM-3; SPTCL; TIMD3; Tim-3; TIMD-3; HAVcr-2
Species	Human
Gene ID (Human)	84868
Gene ID	84868
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-200 of human TIM-3/HAVCR2 (NP_116171.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8TDQ0
Protein Weight	33kDa
Shipping	Ice bag