


Rabbit Anti-Human Tlr1 Monoclonal Antibody For Fc - Q15399

Item Number: CM0010184



Introduction

-conjugated rabbit monoclonal antibody targeting human TLR1 (CD281), validated for flow cytometry. Recognizes toll-like receptor 1 involved in innate immune response to bacterial lipopeptides. Suitable for immunology and infectious disease research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human TLR1 mAb
Remarks/Alias	TIL; CD281; rsc786; TIL. LPRS5
Species	Human
GeneID (Human)	7096
GeneID	7096
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 325-424 of human TLR1 (NP_003254.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q15399
Protein Weight	90kDa
Shipping	Ice bag