


Rabbit Anti-Human Tlr1 Monoclonal Antibody For Fc - Q15399

Item Number: CM0008446



Introduction

-conjugated rabbit monoclonal antibody against human Toll-like receptor 1 (TLR1/CD281). Validated for flow cytometry (FC), suitable for innate immunity research on TLR2/TLR1 heterodimer-mediated lipopeptide recognition.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human TLR1 mAb
Remarks/Alias	TIL; CD281; rsc786; TIL. LPRS5
Species	Human
GeneID (Human)	7096
GeneID	7096
Immunogen	Synthetic peptide (a synthetic peptide corresponding to a sequence within amino acids 325-424 of human TLR1 (NP_003254.2))
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q15399
Protein Weight	90kDa
Shipping	Ice bag