

Rabbit Anti-Human Tnf-A Monoclonal Antibody For Flow Cytometry (Intracellular) - P01375

Item Number: CM0005029

Introduction

-conjugated rabbit monoclonal antibody targeting human TNF- α , validated for intracellular flow cytometry. Specific for human TNF- α , 26 kDa protein. Ideal for immuno-oncology and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human TNF- α mAb
Remarks/Alias	DIF; TNFA; TNFSF2; TNLG1F; TNF-alpha
Species	Human
Gene ID (Human)	7124
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 77-233 of human TNF- α (NP_002344.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P01375
Protein Weight	26kDa
Shipping	Ice bag