


Rabbit Anti-Human Tnf-A Monoclonal Antibody For Intracellular Flow Cytometry - P01375

Item Number: CM0004268



Introduction

-conjugated rabbit monoclonal antibody against human TNF- α , validated for intracellular flow cytometry. Detects a 26 kDa protein. Ideal for cytokine research and immune response studies.

[Learn More](#)

Parameter	Details
Product Name	Rabbit anti-Human TNF- α mAb
Remarks/Alias	DIF; TNFA; TNFSF2; TNLG1F; TNF-alpha
Species	Human
Gene ID (Human)	7124
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 77-233 of human TNF- α (NP_002344.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P01375
Protein Weight	26kDa
Shipping	Ice bag