

Rabbit Anti-Human Tnfrsf10D/Trailr4 Monoclonal Antibody For Fc - Q9Ubn6

Item Number: CM0012262

Introduction

-conjugated rabbit monoclonal antibody targeting human TNFRSF10D/TRAILR4 (CD264, DCR2), a decoy receptor that blocks TRAIL-mediated apoptosis. Suitable for flow cytometry (FC) and validated on human cells. Ideal for apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human TNFRSF10D/TRAILR4 mAb
Remarks/Alias	DCR2; CD264; TRUNDD; TRAILR4; TRAIL-R4
Species	Human
GeneID (human)	8793
GeneID	8793
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 56-211 of human TNFRSF10D/TRAILR4 (NP_003831.2))
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9UBN6
Protein Weight	42kDa
Shipping	Ice bag