

Rabbit Anti-Human Trem1 Monoclonal Antibody For Fc - Q9Np99

Item Number: CM0009396



Introduction

-conjugated rabbit monoclonal antibody specific for human TREM1 (CD354) validated for flow cytometry. Suitable for studying innate immunity, inflammatory responses, and myeloid cell signaling.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human TREM1 mAb
Remarks/Alias	CD354; TREM-1
Species	Human
Gene ID (Human)	54210
Gene ID	54210
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-200 of human TREM1 (NP_061113.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9NP99
Protein Weight	18kDa/26kDa
Shipping	Ice bag