


Rabbit Anti-Human Trop-2 Monoclonal Antibody () For Flow Cytometry - P09758

Item Number: CM0009588



Introduction

Rabbit monoclonal anti-human TROP-2 antibody conjugated to . Detects tumor-associated calcium signal transducer 2 by flow cytometry. Expressed in placenta and pancreatic carcinomas. Ideal for IVD research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human TROP-2 mAb
Remarks/Alias	EGP1; GP50; M1S1; EGP-1; TROP2; GA7331; GA733-1
Species	Human
GeneID	4070
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-274 of human (NP_002344.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P09758
Protein Weight	35kDa
Shipping	Ice bag