

Rabbit Anti-Human Trop-2 Monoclonal Antibody For Flow Cytometry - P09758

Item Number: CM0010595



Introduction

-conjugated rabbit monoclonal antibody targeting human TROP-2 (TACSTD2) for flow cytometry. Ideal for IVD assay development and tumor antigen research.

[Learn More](#)

Parameter	Details
Product Name	Rabbit anti-Human TROP-2 mAb
Remarks/Alias	TACSTD2; EGP-1; EGP1; GA733-1; GA7331; GP50; M1S1; TROP2; tumor-associated calcium signal transducer 2
Species	Human
Gene ID (Human)	4070
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-274 of human TROP-2 (NP_002344.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P09758
Protein Weight	35kDa
Shipping	Ice bag