

Rabbit Anti-Human Trop-2 Monoclonal Antibody For If/Icc, Fc - P09758

Item Number: CM0009134



Introduction

Rabbit monoclonal antibody against human TROP-2 (TACSTD2), conjugated to . Validated for immunofluorescence and flow cytometry. Ideal for tumor marker research and diagnostic assay development.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human TROP-2 mAb
Remarks/Alias	EGP1; GP50; M1S1; EGP-1; TROP2; GA7331; GA733-1
Species	Human
Gene ID (Human)	4070
Gene ID	4070
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 27-274 of human TROP-2 (NP_002344.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Human
SWISS	P09758
Protein Weight	36kDa
Shipping	Ice bag