

Rabbit Anti-Human Ulbp2 Monoclonal Antibody For Fc - Q9Bzm5

Item Number: CM0011443



Introduction

-conjugated rabbit monoclonal antibody targeting human ULBP2 (NKG2DL2). Designed for flow cytometry applications in cancer and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human ULBP2 mAb
Remarks/Alias	N2DL2; RAET1H; RAET1L; NKG2DL2; ALCAN-alpha
Species	Human
GeneID (human)	80328
GeneID	80328
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-217 of human ULBP2 (NP_079493.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9BZM5
Protein Weight	27kDa
Shipping	Ice bag