

Rabbit Anti-Human B2 Microglobulin Mab For Flow Cytometry - P61769

Item Number: CM0009481



Introduction

-conjugated rabbit anti-human β 2-microglobulin monoclonal antibody for flow cytometry detection of B2M in immune monitoring and MHC expression studies. This reagent targets the invariant chain of MHC class I, enabling the study of antigen presentation and immune function.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human β 2 Microglobulin mAb
Remarks/Alias	IMD43
Species	Human
GeneID (Human)	567
GeneID	567
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-119 of human β 2 Microglobulin (NP_004039.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P61769
Protein Weight	14kDa
Shipping	Ice bag