

Rabbit Anti-Human B2 Microglobulin Mab - P61769

Item Number: CM0012369



Introduction

-conjugated rabbit monoclonal antibody targeting human β 2 Microglobulin (B2M), recommended for flow cytometry. Suitable for immune system research, MHC class I studies, and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human β 2 Microglobulin mAb
Remarks/Alias	B2M; IMD43; beta-2-microglobulin
Species	Human
GeneID	567
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-119 of human β 2 Microglobulin (NP_004039.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P61769
Protein Weight	13kDa
Shipping	Ice bag