

Rabbit Anti-Human B2 Microglobulin Monoclonal Antibody For Flow Cytometry - P61769

Item Number: CM0009298



Introduction

Rabbit monoclonal anti-human β 2 Microglobulin antibody (clone , IMD43) validated for flow cytometry. Specific to human; recognizes the 14kDa B2M protein. Ideal for immune research and IVD applications.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human β 2 Microglobulin mAb
Remarks/Alias	IMD43
Species	Human
Gene ID	567
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 21-119 of human β 2 Microglobulin (NP_004039.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P61769
Protein Weight	14kDa
Shipping	Ice bag