

Rabbit Anti-Human B2 Microglobulin Monoclonal Antibody For Flow Cytometry - P61769

Item Number: CM0011428



Introduction

-conjugated rabbit monoclonal antibody targeting human β 2-Microglobulin (B2M). Validated for flow cytometry (FC). Supports immune research, MHC class I antigen presentation studies, and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human β 2 Microglobulin mAb
Remarks/Alias	B2M; IMD43; beta-2-microglobulin
Species	Human
GeneID (Human)	567
GeneID	567
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-119 of human β 2 Microglobulin (NP_004039.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P61769
Protein Weight	13kDa
Shipping	Ice bag