

Rabbit Anti-Human B2 Microglobulin Monoclonal Antibody For Flow Cytometry - P61769

Item Number: CM0009230



Introduction

-conjugated rabbit monoclonal antibody targeting human $\beta 2$ Microglobulin for flow cytometry. Detects a 14kDa secreted and cell surface protein involved in MHC class I antigen presentation. Suitable for immunological research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human $\beta 2$ Microglobulin mAb
Remarks/Alias	IMD43
Species	Human
GeneID	567
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-119 of human $\beta 2$ Microglobulin (NP_004039.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P61769
Protein Weight	14kDa
Shipping	Ice bag