


Rabbit Anti-Human Vcam-1/Cd106 Mab For Fc - P19320

Item Number: CM0004837



Introduction

Rabbit anti-Human VCAM-1/CD106 monoclonal antibody for flow cytometry (FC). Validated for human samples; targets vascular cell adhesion protein 1 involved in immune surveillance and inflammation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human VCAM-1/CD106 mAb
Remarks/Alias	CD106; INCAM-100
Species	Human
Gene ID	7412
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-698 of human VCAM-1/CD106 (NP_001069.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P19320
Protein Weight	71kDa/74kDa/81kDa
Shipping	Ice bag