

Rabbit Anti-Human Vcam-1/Cd106 Monoclonal Antibody (Fc) - P19320

Item Number: CM0006092



Introduction

-conjugated rabbit monoclonal antibody targeting human VCAM-1 (CD106). Validated for flow cytometry (FC). Recombinant immunogen (aa 25-698). Ideal for inflammation and leukocyte adhesion research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human VCAM-1/CD106 mAb
Remarks/Alias	CD106; INCAM-100
Species	Human
Gene ID (Human)	7412
Gene ID	7412
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 25-698 of human VCAM1 (NP_001069.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P19320
Protein Weight	71kDa/74kDa/81kDa
Shipping	Ice bag