

Rabbit Anti-Human Vcam-1/Cd106 Monoclonal Antibody For Fc - P19320

Item Number: CM0005352



Introduction

Recombinant rabbit monoclonal antibody targeting human VCAM-1 (CD106; INCAM-100). Validated for flow cytometry (FC). Detects multiple isoforms (~71, 74, 81 kDa). Suitable for research on leukocyte adhesion, inflammation, and immune surveillance.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human VCAM-1/CD106 mAb
Remarks/Alias	CD106; INCAM-100
Species	Human
GeneID (Human)	7412
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-698 of human VCAM-1/CD106 (NP_001069.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P19320
Protein Weight	71kDa/74kDa/81kDa
Shipping	Ice bag