

Rabbit Anti-Human/Monkey Cd127/IL-7R α Mab - P16871

Item Number: CM0010828



Introduction

-conjugated rabbit monoclonal antibody against human and cynomolgus CD127/IL-7R α for flow cytometry. Ideal for immunology research and diagnostic assay development involving IL-7 receptor pathways.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD127/IL-7R α mAb
Remarks/Alias	ILRA; CD127; IL7RA; CDW127; IMD104; sIL-7R; Inc-IL7R; IL7Ralpha; IL-7Ralpha; IL-7R-alpha; CD127/IL7R
Species	Human
GeneID (human)	3575
GeneID	3575
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P16871
Protein Weight	52kDa
Shipping	Ice bag