

Rabbit Anti-Human/Monkey Cd127/IL-7R α Monoclonal Antibody For Flow Cytometry - P16871

Item Number: CM0012005



Introduction

Rabbit anti-Human/Monkey CD127/IL-7R α Monoclonal Antibody validated for Flow Cytometry. Recombinant immunogen, reacts with Human and Cynomolgus. Ideal for IL-7R expression analysis in immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD127/IL-7R α mAb
Remarks/Alias	ILRA; CD127; IL7RA; CDW127; IMD104; sIL-7R; Inc-IL7R; IL7Ralpha; IL-7Ralpha; IL-7R-alpha; CD127/IL7R
Species	Human
GeneID (Human)	3575
GeneID	3575
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-236 human CD127/IL-7R α (NP_002176.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P16871
Protein Weight	52kDa
Shipping	Ice bag