


Rabbit Anti-Human/Monkey Cd14 Mab For Fc - P08571

Item Number: CM0010761



Introduction

-conjugated rabbit monoclonal antibody targeting human and cynomolgus CD14 for flow cytometry. CD14 is a GPI-anchored LPS coreceptor mediating innate immunity. Ideal for monocyte/macrophage identification in diagnostic research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD14 mAb
Remarks/Alias	CD14
Species	Human
Gene ID (Human)	929
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-352 of human CD14 (NP_000582.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P08571
Protein Weight	40kDa
Shipping	Ice bag