


Rabbit Anti-Human/Monkey Cd14 Monoclonal Antibody For Fc - P08571

Item Number: CM0008087



Introduction

-conjugated rabbit monoclonal antibody specific to human and cynomolgus CD14, validated for flow cytometry. Targets monocyte differentiation antigen CD14, a coreceptor for bacterial lipopolysaccharide and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD14 mAb
Remarks/Alias	CD14
Species	Human
GeneID (Human)	929
GeneID	929
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 20-352 of human CD14 (NP_000582.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P08571
Protein Weight	40kDa
Shipping	Ice bag