


Rabbit Anti-Human/Monkey Cd14 Monoclonal Antibody For Flow Cytometry - P08571

Item Number: CM0012520



Introduction

-conjugated rabbit monoclonal antibody targeting CD14 (human/cynomolgus). Designed for flow cytometry (FC) to identify monocytes/macrophages and study innate immune responses. P08571.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD14 mAb
Remarks/Alias	CD14
Species	Human
GeneID	929
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-352 of human CD14 (NP_000582.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P08571
Protein Weight	40kDa
Shipping	Ice bag