

Rabbit Anti-Human/Monkey Cd14 Monoclonal Antibody For Flow Cytometry - P08571

Item Number: CM0007878



Introduction

Rabbit anti-Human/Monkey CD14 Monoclonal Antibody validated for Flow Cytometry. Recognizes CD14, a LPS coreceptor mediating innate immune responses. Suitable for human and cynomolgus samples. Ideal for immunology and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD14 mAb
Remarks/Alias	CD14
Species	Human
GeneID	929
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-352 of human CD14 (NP_000582.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P08571
Protein Weight	40kDa
Shipping	Ice bag