

Rabbit Anti-Human/Monkey Cd16 Mab For Fc - P08637/O75015

Item Number: CM0008981



Introduction

Rabbit anti-Human/Monkey CD16 monoclonal antibody for flow cytometry. Recognizes CD16A and CD16B (FCGR3A/FCGR3B), with reactivity in human and cynomolgus monkey. Suitable for FC application.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD16 mAb
Remarks/Alias	CD16; FCG3; CD16A; FCGR3; IGFR3; IMD20; FCR-10; FCRIII; CD16-II; FCGRIII; FCRIIIA; FcGRIIIA; CD16b
Species	Human
Gene ID (Human)	2214/2215
Immunogen	Recombinant protein/Recombinant fusion protein containing a sequence corresponding to amino acids 20-208 of human FCGR3A/FCGR3B (P08637/O75015)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P08637/O75015
Protein Weight	29KDa/26KDa
Shipping	Ice bag