

Rabbit Anti-Human/Monkey Cd16 Monoclonal Antibody For Flow Cytometry - P08637/O75015

Item Number: CM0006131



Introduction

-conjugated rabbit monoclonal anti-CD16 antibody validated for flow cytometry, targeting human and cynomolgus CD16 (FCGR3A/FCGR3B). Recombinant immunogen covers extracellular domain, ensuring specific detection of low-affinity IgG Fc receptor involved in ADCC.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD16 mAb
Remarks/Alias	CD16; FCG3; CD16A; FCGR3; IGFR3; IMD20; FCR-10; FCRIII; CD16-II; FCGRIII; FCRIIIA; FcGRIIIA; CD16b
Species	Human
Gene ID (Human)	2214/2215
Gene ID	2214/2215
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-208 of human FCGR3A/FCGR3B (P08637/O75015).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P08637/O75015
Protein Weight	29KDa/26KDa
Shipping	Ice bag