

Rabbit Anti-Human/Monkey Cd18 Mab - P05107

Item Number: CM0004059



Introduction

-conjugated rabbit monoclonal antibody targeting human and cynomolgus CD18 for flow cytometry. CD18 (ITGB2) is the beta-2 integrin subunit mediating leukocyte adhesion and migration.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD18 mAb
Remarks/Alias	LAD; CD18; MF17; MFI7; LCAMB; LFA-1; MAC-1
Species	human
GeneID (Human)	3689
GeneID	3689
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-769 of human CD18 (NP_000202.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P05107
Protein Weight	85kDa
Shipping	Ice bag