

Rabbit Anti-Human/Monkey Cd19 Mab - Cd19

Item Number: CM0012338

Introduction

Anti-CD19 monoclonal antibody (clone) produced in rabbit, validated for flow cytometry. Recognizes human and cynomolgus CD19, a 95 kDa B-cell surface glycoprotein. Suitable for hematopoietic and immunological research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD19 mAb
Remarks/Alias	CD19
Species	Monkey
Gene ID (Human)	930
Gene ID	102145514
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-292 of monkey CD19 (XP_005591597.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
Shipping	Ice bag