

Rabbit Anti-Human/Monkey Cd19 Mab - Cd19

Item Number: CM0011139

Introduction

Rabbit anti-Human/Monkey CD19 monoclonal antibody for flow cytometry (FC). Cross-reacts with human and cynomolgus monkey CD19 (63kDa). Ideal for B-cell research and immunophenotyping.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD19 mAb
Alias	CD19
Species	Monkey
Gene ID (Human)	930
Gene ID (Monkey)	102145514
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-292 of monkey CD19
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
Protein Weight	63kDa
Shipping	Ice bag