

Rabbit Anti-Human/Monkey Cd19 Mab - Cd19

Item Number: CM0009318

Introduction

-conjugated rabbit monoclonal antibody targeting CD19, reactive with human and cynomolgus monkey. Designed for flow cytometry applications, this antibody recognizes a 63 kDa transmembrane protein, making it suitable for B-cell immunophenotyping research.

[Learn More](#)

Specification	Details
Product Name	Rabbit anti-Human/Monkey CD19 mAb
Remarks/Alias	CD19
Species	Monkey
Gene ID (Human)	930
Gene ID	102145514
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-292 of monkey CD19 (XP_005591597.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
Protein Weight	63kDa
Shipping	Ice bag