

Rabbit Anti-Human/Monkey Cd19 Mab

Item Number: CM0005132

Introduction

-conjugated rabbit monoclonal antibody targeting human and monkey CD19. Validated for flow cytometry with cross-reactivity to human and cynomolgus monkey. Predicted molecular weight: 63 kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD19 mAb
Species	Monkey
Gene ID (Human)	930
Gene ID	102145514
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-292 of monkey CD19
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
Protein Weight	63kDa
Shipping	Ice bag