

Rabbit Anti-Human/Monkey Cd19 Monoclonal Antibody - Cd19

Item Number: CM0011307

Introduction

-conjugated rabbit monoclonal antibody targeting CD19, validated for flow cytometry. Reacts with human and cynomolgus monkey CD19, aiding B cell research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD19 mAb
Remarks/Alias	CD19
Species	Monkey
Gene ID (Human)	930
Gene ID	102145514
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-292 of monkey CD19
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
Shipping	Ice bag