

Rabbit Anti-Human/Monkey Cd19 Monoclonal Antibody - Cd19

Item Number: CM0003981

Introduction

Recombinant -conjugated rabbit monoclonal antibody targeting CD19 for flow cytometry. Validated for human and cynomolgus monkey species. Suitable for B-cell lineage analysis.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD19 mAb
Remarks/Alias	CD19
Species	Monkey
Gene ID (Human)	930
Gene ID (Monkey)	102145514
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-292 of monkey CD19.
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
Protein Weight	63kDa
Shipping	Ice bag