

## Rabbit Anti-Human/Monkey Cd20 Mab For Fc - P12314

Item Number: CM0004676

### Introduction

High affinity rabbit monoclonal antibody against human/cynomolgus CD20, conjugated with , validated for flow cytometry. Ideal for immune cell profiling and clinical research.

[Learn More](#)

| Parameter        | Value                                   |
|------------------|-----------------------------------------|
| Product Name     | Rabbit anti-Human/Monkey CD20 mAb       |
| Remarks/Alias    | B1; S7; Bp35; CD20; FMC7; CVID5; LEU-16 |
| Species          | Human                                   |
| Gene ID (Human)  | 931                                     |
| Gene ID          | 931                                     |
| Immunogen        | Synthetic peptide                       |
| Source           | Rabbit                                  |
| Category         | Monoclonal Antibodies                   |
| Application      | FC                                      |
| Cross Reactivity | Human, Cynomolgus                       |
| SWISS            | P12314                                  |
| Protein Weight   | 14kDa/33kDa                             |
| Shipping         | Ice bag                                 |