


Rabbit Anti-Human/Monkey Cd20 Mab For Flow Cytometry - P11836

Item Number: CM0005246



Introduction

-conjugated rabbit monoclonal antibody targeting human/monkey CD20, optimized for flow cytometry. Recognizes B-cell specific 33/14 kDa isoforms; cross-reacts with cynomolgus monkey.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human/Monkey CD20 mAb
Remarks/Alias	B1; S7; Bp35; CD20; FMC7; CVID5; LEU-16
Species	Human
Gene ID (Human)	931
Gene ID	931
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P11836
Protein Weight	14kDa/33kDa
Shipping	Ice bag